



## Overview

- Rotating padded seat
- Manual adjustable base
- Braked rear castors
- Excellent substitution for a wheelchair
- Maximum dimensions: 105 cm (height), 93 cm (length), 64 cm (width)
- Manual adjustment of base width (three levels):
  - Between castors: 45 cm (narrow), 64 cm (middle), 83 cm (wide)
  - Outside of castors: 58 cm (narrow), 77 cm (middle), 97 cm (wide)
- Castors with diameter of 75 mm, rear
- Castors are equipped with brake
- Seamless padding 5 cm
- Trolley weight: 33 kg
- Safe working load: 150 kg