



# Service & Repair





## The Team

**Our nationwide service and Installation teams specialise in onsite LOLER, PUWER and PAT testing across all moving, handling and pressure area care equipment.**

All our employees are Dementia Friends and have been specifically trained by manufacturers to install, test and repair their equipment. This helps to realise a strong return on investment for each asset, increasing the in-life use by up to 57%.

We adhere to the strict guidelines set out by the CQC, HSE and MHRA to ensure your organisation complies with the required standards and safety to protect you, your team and the residents.



## Service & Repair

**Delivered Health Solutions can ensure all your equipment adheres to:**

**LOLER** (Lifting Operations and Lifting Equipment Regulations 1998)

**PUWER** (Provision and Use Of Work Equipment Regulations 1998)

**PAT Testing** (Portable Appliance Test)

Our team of highly skilled engineers will perform your 6-monthly LOLER & 12 monthly PUWER tests across all the relevant moving and handling equipment:

If your equipment fails the test, Delivered Health Solutions will endeavour to repair it on-site with our comprehensive range of spares on our vehicles. If this is not possible we can often recycle parts from other condemned machines which keeps the costs down.

Where this is not possible, we will instruct our workshop to send out the relevant parts and return to fix the issue as soon as possible.

Where parts are no longer available or the product is damaged beyond repair, we have strong partnerships across all manufacturers to find you a replacement at the right price.



LOLER and PUWER Servicing & Repair  
Installation and Training | Asset Management

Care Shop | Unit B Circle Court | Trafford Park | Warren Bruce Road | Manchester, M17 1LB

[careshop.co.uk](http://careshop.co.uk)

0161 848 4000

[sales@deliveredhealthsolutions.co.uk](mailto:sales@deliveredhealthsolutions.co.uk)